

# RASSF7 antibody



Catalog Number: 114484

## Product name

RASSF7 antibody

## Immunogen

[Human RASSF7 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RASSF7 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 36kd,40kd

## Preparation

This antibody was obtained by immunization of RASSF7 recombinant protein (Accession Number: XM\_017018361). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114484(RASSF7 antibody) at dilution of 1:200