

Anti-Nuclear Matrix Protein p84 Rabbit antibody



Catalog Number: 167484

Product name

WB, ICC/IF

Anti-Nuclear Matrix Protein p84 Rabbit antibody

Dilutions

Specificity

WB: 1/1000-1/5000

Human

ICC/IF: 1/20-1/50

Antibody description

Nuclear Matrix Protein p84 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Nuclear Matrix Protein p84

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Western blot detection of Nuclear Matrix Protein p84 in Jurkat, HeLa, COS7 cell lysates using Nuclear Matrix Protein p84 Rabbit pAb(1:1000 diluted). Predicted band size: 76KDa. Observed band size: 84KDa.



Immunocytochemistry of Nuclear Matrix Protein p84(green) in HeLa cells using Nuclear Matrix Protein p84 Rabbit pAb at dilution 1/50, and DAPI(blue)