

YES1-Specific antibody



Catalog Number: 116884

Product name

YES1-Specific antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

YES1-Specific Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue, HeLa cells. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 50-60 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005433). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116884(YES1-Specific antibody) at dilution of



Immunohistochemical of paraffin-embedded human colon using Catalog No:116884(YES1-Specific antibody) at dilution of 1:50 (under 10x



Immunohistochemical of paraffin-embedded


YES1-Specific antibody



Catalog Number: 116884

human colon using Catalog No:116884(YES1-Specific antibody) at dilution of 1:50 (under 40x

IP Result of anti-YES1-Specific (IP:Catalog No:116884, 4ug; Detection:Catalog No:116884 1:300) with human placenta tissue lysate 1520ug.



lens)