

# Anti-Fetuin-A antibody



Catalog Number: 106184

## Product name

Anti-Fetuin-A antibody

## Immunogen

[Mouse Fetuin-A \(His Tag\) recombinant protein](#)

## Specificity

Mouse Fetuin-A / AHSG

**Has cross-reactivity** in ELISA with ; Human FetuinA/AHSG

## Antibody description

Rabbit monoclonal to Fetuin-A

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Fetuin-A / AHSG (rM Fetuin-A / AHSG; NP\_032037.2; Met1-Ile345).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



Fetuin-A / AHSG Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse AHSG in mouse liver (from 2 donors) with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).