

Anti-EGFR Rabbit antibody



Catalog Number: 167184

Product name

WB, ICC/IF, FC

Anti-EGFR Rabbit antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse

ICC/IF: 1/50

Antibody description

FC: 1/20

EGFR Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human EGFR

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of EGFR(green) in HeLa cells using EGFR Rabbit pAb at dilution 1/50, and



DAPI(blue)

Western blot detection of EGFR in HeLa cell lysates using EGFR Rabbit pAb(1:1000 diluted). Predicted band size: 134kDa. Observed band size: 175kDa.