

CXCL16 antibody



Catalog Number: 107184

Product name

CXCL16 antibody

Clonality

Monoclonal

Immunogen

[Human CXCL16 Recombinant protein \(His tag\)](#)

Ig Type

Mouse IgG2a

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

CXCL16 Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 35 kDa, 60 kDa

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of CXCL16 recombinant protein (Accession Number: NM_022059). Purification method: Caprylic acid/ammonium sulfate precipitation.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107184(CXCL16 antibody) at dilution of 1:300