

Product name

SCAP antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse; other species not tested.

Antibody description

SCAP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, K-562 cells. Positive IHC detected in human lung cancer tissue, human kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 140kd,96kd

Preparation

This antibody was obtained by immunization of SCAP recombinant protein (Accession Number: BC020987). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

Validations

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114984(SCAP antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:114984(SCAP Antibody) at dilution of 1:50 (under 40x lens)




Immunohistochemistry of paraffin-embedded

SCAP antibody



Catalog Number: 114984

human kidney tissue slide using Catalog No:114984(SCAP Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HEK-293 cells using Catalog No:114984(SCAP Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)