

# ALKBH1 antibody



Catalog Number: 107984

## Product name

ALKBH1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ALKBH1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human colon tissue, human pancreas tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of ALKBH1 recombinant protein (Accession Number: NM\_006020). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107984(ALKBH1 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107984(ALKBH1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107984(ALKBH1 Antibody) at dilution of 1:50

# ALKBH1 antibody



Catalog Number: 107984

---

(under 40x lens)