

NSBP1 antibody



Catalog Number: 113284

Product name

NSBP1 antibody

Specificity

Human; other species not tested.

Antibody description

NSBP1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in PC-3 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of NSBP1 recombinant protein (Accession Number: NM_030763). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113284(NSBP1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113284(NSBP1 Antibody) at dilution of 1:200 (under 40x lens)