

# PCBD2 antibody



Catalog Number: 113684

## Product name

PCBD2 antibody

## Immunogen

[Human PCBD2 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PCBD2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 14 kDa

## Preparation

This antibody was obtained by immunization of PCBD2 recombinant protein (Accession Number: NM\_032151). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113684(PCBD2 antibody) at dilution of 1:500