

Ephrin A3 antibody



Catalog Number: 110384

Product name

Ephrin A3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Ephrin A3 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A549 cells. Positive IP detected in A549 cells. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 26-35 kDa

Preparation

This antibody was obtained by immunization of Ephrin A3 recombinant protein (Accession Number: NM_004952). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110384(EFNA3 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:110384(EFNA3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded


Ephrin A3 antibody



Catalog Number: 110384

human skeletal muscle using Catalog
No:110384(EFNA3 antibody) at dilution of 1:50

IP Result of anti-EFNA3 (IP:Catalog No:110384,
4ug; Detection:Catalog No:110384 1:500) with
A549 cells lysate 1040ug.



(under 40x lens)