

TTC39A antibody



Catalog Number: 116384

Product name

TTC39A antibody

Immunogen

[Human TTC39A Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

TTC39A Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, MCF7 cells. Positive IF detected in HepG2. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 65-67kd

Preparation

This antibody was obtained by immunization of TTC39A recombinant protein (Accession Number: NM_001297666). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116384(TTC39A antibody) at dilution of 1:100




Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116384(TTC39A antibody) at dilution of 1:50 (under 10x lens)

TTC39A antibody

Catalog Number: 116384



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116384(TTC39A antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using TTC39A antibody Catalog No:116384 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).