

Anti-Alpha Skeletal Muscle Actin Rabbit antibody



Catalog Number: 166784

Product name

Anti-Alpha Skeletal Muscle Actin Rabbit antibody

IP: 1/20

Specificity

Human, Mouse, Rat

Antibody description

Alpha Skeletal Muscle Actin Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Actin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions


WB: 1/1000

IHC: 1/50


ICC/IF: 1/20

FC: 1/20


Validations



Western blot detection of Actin in Jurkat, Ramos, MCF7 cell lysates using Actin Rabbit pAb(1:1000 diluted). Predicted band size: 42KDa. Observed band size: 42KDa.



Immunohistochemistry of Actin in paraffin-embedded Human colon cancer tissue using Actin



Rabbit pAb at dilution 1/100

Western blot detection of Actin in HeLa cell lysates using Actin Rabbit pAb(1:1000 diluted). Predicted band size: 42KDa. Observed

Anti-Alpha Skeletal Muscle Actin Rabbit antibody



Catalog Number: 166784

band size:42KDa.