

PAIP1 antibody



Catalog Number: 113494

Product name

PAIP1 antibody

Immunogen

[Human PAIP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat, Monkey; other species not tested.

Antibody description

PAIP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, BxPC-3 cells, COS-7 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 60 kDa, 51 kDa, 40 kDa

Preparation

This antibody was obtained by immunization of PAIP1 recombinant protein (Accession Number: NM_182789). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113494(PAIP1 antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using PAIP1 antibody Catalog No:113494 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).