

MUL1 antibody



Catalog Number: 112894

Product name

MUL1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MUL1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human spleen tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of MUL1 recombinant protein (Accession Number: NM_024544). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112894(MUL1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112894(MUL1 antibody) at dilution of 1:50 (under 40x lens)




Immunofluorescent analysis of HepG2 cells, using Catalog No:112894 and Rhodamine-labeled goat

MUL1 antibody



Catalog Number: 112894

anti-rabbit IgG (red).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112894(MUL1 Antibody) at dilution of 1:1000