

Anti-Phospho-EGFR (Tyr1173) Rabbit antibody



Catalog Number: 166994

Product name

Anti-Phospho-EGFR (Tyr1173) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Phospho-EGFR (Tyr1173) Rabbit polyclonal antibody

Dilutions

WB: 1/500

Preparation

Antigen: CSLDNPDP-YQQDFFP

Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Validations



Western blot analysis of extracts of A549 cells, untreated or EGF-stimulated, using Phospho-EGF Receptor (Tyr1173) Rabbit pAb (upper, 1:500 diluted) and EGF Receptor Antibody 201012-6H11 (lower, 1:500 diluted).