

# Anti-Top2a antibody



Catalog Number: 176694

## Product name

WB, ICC

Anti-Top2a antibody

## Dilutions

## Specificity

WB: 1:1,000-1:2,000

Human, Mouse, Rat

ICC: 1:200

## Antibody description

Rabbit polyclonal antibody to Top2a

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

## Validations




Fig1: Western blot analysis on Jurkat cell lysates using anti-TOP2A rabbit polyclonal antibodies.




Fig2: ICC staining TOP2A in A431 cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS and counterstained with DAPI in order to highlight the nucleus (blue).