

ACSF2 antibody



Catalog Number: 107694

Product name

ACSF2 antibody

Immunogen

[Human ACSF2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACSF2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IHC detected in human skin tissue. Observed molecular weight by Western-blot: 55-70kd

Preparation

This antibody was obtained by immunization of ACSF2 recombinant protein (Accession Number: NM_025149). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107694(ACSF2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:107694(ACSF2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:107694(ACSF2 Antibody) at dilution of 1:50 (under 40x lens)