

TFCP2 antibody



Catalog Number: 116005

Product name

TFCP2 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

TFCP2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 66kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of TFCP2 recombinant protein (Accession Number: NM_005653). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116005(TFCP2 antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality