

ZNF571 antibody



Catalog Number: 117005

Product name

ZNF571 antibody

Recommended Dilution:

WB: 1:200-1:1000

Immunogen

IP: 1:200-1:1000

[Human ZNF571 Recombinant protein \(His tag\)](#)

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

ZNF571 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive IP detected in MCF-7 cells. Positive WB detected in MCF-7 cells, HeLa cells. Observed molecular weight by Western-blot: 65 kDa




Preparation

This antibody was obtained by immunization of ZNF571 recombinant protein (Accession Number: NM_001290314). Purification method: Antigen affinity purified.

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117005(ZNF571 Antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

IP Result of anti-ZNF571 (IP:Catalog No:117005, 4ug; Detection:Catalog No:117005 1:300) with MCF-7 cells lysate 1200ug.

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions

ZNF571 antibody



Catalog Number: 117005



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:117005(ZNF571 Antibody) at dilution of 1:50 (under 40x lens)