

# GABPB1 antibody



Catalog Number: 110805

## Product name

GABPB1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GABPB1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HeLa cells, HepG2 cells, MCF7 cells, mouse testis tissue. Positive IP detected in HeLa cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 36-45 kDa, 47-53 kDa

## Preparation

This antibody was obtained by immunization of GABPB1 recombinant protein (Accession Number: NM\_002041). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, WB, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110805(GABPB1 antibody) at dilution of



Immunohistochemical of paraffin-embedded human testis using Catalog No:110805(GABPB1 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-GABPB1 (IP:Catalog No:110805,

# GABPB1 antibody



Catalog Number: 110805

---

4ug; Detection:Catalog No:110805 1:1000) with HeLa cells lysate 520ug.