

PHD3 antibody



Catalog Number: 113805

Product name

PHD3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

PHD3 Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells, A375 cells, HEK-293 cells, Jurkat cells. Positive IP detected in A375 cells. Positive IHC detected in human brain tissue, human heart tissue. Observed molecular weight by Western-blot: 28-32k kDa

Preparation

This antibody was obtained by immunization of PHD3 recombinant protein (Accession Number: NM_022073). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:113805(EGLN3 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:113805(EGLN3 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:113805(EGLN3

PHD3 antibody



Catalog Number: 113805

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-EGLN3 (IP:Catalog No:113805, 4ug; Detection:Catalog No:113805 1:500) with A375 cells lysate 3600ug.