

RAC2 antibody



Catalog Number: 114505

Product name

RAC2 antibody

Specificity

human, mouse; other species not tested.

Antibody description

RAC2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HL-60 cells. Observed molecular weight by Western-blot: 21kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002872). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114505(RAC2 antibody) at dilution of 1:100