

ALX4 Polyclonal Antibody



Catalog Number: 163905

Product name

Rabbit IgG

ALX4 Polyclonal Antibody

Applications

Specificity

WB, IHC

Human, Mouse, Rat

Dilutions

Antibody description

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Polyclonal antibody to ALX4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ALX4 (NP_068745.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - ALX4 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using ALX4 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.