

Anti-RAB7 Rabbit antibody



Catalog Number: 167205

Product name

WB: 1/1000-1/5000

Anti-RAB7 Rabbit antibody

IHC: 1/20

Specificity

ICC/IF: 1/20

Human, Mouse, Rat

FC: 1/20

Antibody description

RAB7 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human RAB7

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Western blot detection of RAB7 in 3T3 cell lysates using RAB7 Rabbit pAb(1:1000 diluted). Predicted band size:24kDa. Observed band size:24kDa.



Immunohistochemistry of RAB7 in paraffin-embedded Human lung cancer tissue using RAB7 Rabbit pAb at dilution 1/20