

C16orf45 antibody



Catalog Number: 108605

Product name

C16orf45 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

C16orf45 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 24-26 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of C16orf45 recombinant protein (Accession Number: NM_033201). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108605(C16orf45 Antibody) at dilution of 1:600

Clonality