

Product name

NAGA antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NAGA Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 48 kDa, 60-70 kDa

Preparation

This antibody was obtained by immunization of NAGA recombinant protein (Accession Number: NM_001362850). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113015(NAGA antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human placenta using Catalog No:113015(NAGA antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human placenta using Catalog No:113015(NAGA

NAGA antibody



Catalog Number: 113015

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:113015(NAGA Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)