

Product name

MYCBP antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYCBP Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human kidney tissue, human lung tissue, human placenta tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of MYCBP recombinant protein (Accession Number: NM_012333). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112915(MYCBP antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:112915(MYCBP antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of HepG2 cells, using

MYCBP antibody



Catalog Number: 112915

MYCBP antibody Catalog No:112915 at 1:50
dilution and Rhodamine-labeled goat anti-rabbit

IgG (red). Blue pseudocolor = DAPI (fluorescent
DNA dye).