

TAGLN2 Polyclonal Antibody



Catalog Number: 165025

Product name

TAGLN2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TAGLN2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human TAGLN2 (NP_003555.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - TAGLN2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TAGLN2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.


Western blot - TAGLN2 Polyclonal Antibody

Western blot analysis of extracts from normal (control) and TAGLN2 knockout 293T cells, using TAGLN2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 5s.

TAGLN2 Polyclonal Antibody



Catalog Number: 165025



Immunofluorescence - TAGLN2 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using TAGLN2 antibody . Blue: DAPI for nuclear staining.