

CCL25 Polyclonal Antibody



Catalog Number: 164925

Product name

CCL25 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:100 - 1:200

IF 1:50 - 1:100

Specificity

Human

Antibody description

Polyclonal antibody to CCL25

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-150 of human CCL25 (NP_005615.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Immunohistochemistry - CCL25 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using CCL25 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - CCL25 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CCL25 antibody . Blue: DAPI for nuclear staining.