

# Anti-phospho-Cyclin B1 (Ser126) Rabbit antibody



Catalog Number: 167625

## Product name

Monoclonal

Anti-phospho-Cyclin B1 (Ser126) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IP

## Antibody description

phospho-Cyclin B1 (Ser126) Rabbit monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

IP: 1/20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser126 of human Cyclin B1

## Validations


## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality



Western blot detection of Cyclin B1 (phospho-Ser126) in HeLa,Ramos using Cyclin B1 (phospho-Ser126) Rabbit mAb(1:1000 diluted)