

# H3F3B antibody



Catalog Number: 111325

## Product name

H3F3B antibody

## Immunogen

[Human H3F3B Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

H3F3B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 17 kDa

## Preparation

This antibody was obtained by immunization of H3F3B recombinant protein (Accession Number: NM\_001379045). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111325(H3F3B antibody) at dilution of 1:300



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111325(H3F3B antibody) at dilution of 1:300