

# VCAM-1 antibody



Catalog Number: 116725

## Product name

ELISA, IHC

VCAM-1 antibody

## Dilutions

## Specificity

Recommended Dilution:


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

VCAM-1 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human kidney tissue, human lung cancer tissue, human placenta tissue, human tonsillitis tissue.

## Validations




## Preparation

This antibody was obtained by immunization of VCAM-1 recombinant protein (Accession Number: NM\_001078). Purification method: Antigen affinity purified.

Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:116725(VCAM-1 Antibody) at dilution of 1:50 (under 10x lens)

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:116725(VCAM-1 Antibody) at dilution of 1:50 (under 40x lens)

## Applications