

FAAH antibody



Catalog Number: 110435

Product name

FAAH antibody

Immunogen

[Human FAAH Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAAH Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A431 cells, rat testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of FAAH recombinant protein (Accession Number: NM_001441). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110435(FAAH Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110435(FAAH Antibody) at dilution of 1:50 (under 40x lens)