

TRAF7 antibody



Catalog Number: 116235

Product name

TRAF7 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRAF7 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse heart tissue, rat heart tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 67-75 kDa

Preparation

This antibody was obtained by immunization of TRAF7 recombinant protein (Accession Number: BC024267). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations




Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116235(TRAF7 antibody) at dilution of 1:50



(under 10x lens)

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116235(TRAF7 Antibody) at dilution of 1:300



IP Result of anti-TRAF7 (IP:Catalog No:116235, 3ug; Detection:Catalog No:116235 1:300) with

TRAF7 antibody



Catalog Number: 116235

mouse liver tissue lysate 3000ug.