

Anti-UFC1 Rabbit antibody



Catalog Number: 168635

Product name

WB, ICC/IF, IP

Anti-UFC1 Rabbit antibody

Dilutions

Specificity

WB: 1/1000-1/5000

Human, Mouse, Rat

ICC/IF: 1/20

Antibody description

IP: 1/20

UFC1 Rabbit monoclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human UFC1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.


Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications



Immunocytochemistry of UFC1(green) in HeLa cells using UFC1 Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of UFC1 in K562, HeLa cell lysates using UFC1 Rabbit mAb(1:1000 diluted). Predicted band size: 20kDa. Observed band size: 20kDa.