

# Anti-TBPL2 antibody



Catalog Number: 176335

## Product name

Polyclonal

Anti-TBPL2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to TBPL2

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tbpl2. aa range:81-130

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1:: Western Blot analysis of HeLa cells using TBPL2 Polyclonal Antibody