

Anti-ATP6V1A Rabbit antibody



Catalog Number: 167735

Product name

Anti-ATP6V1A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATP6V1A Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human ATP6V1A

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/100

IP: 1/20

Validations



Western blot detection of ATP6V1A in HeLa,A549,HL-60,U251,U87-MG cell lysates using ATP6V1A Rabbit mAb(1:1000 diluted).Predicted band size:68kDa.Observed band size:68kDa.