

CD127/IL7R antibody



Catalog Number: 109045

Product name

CD127/IL7R antibody

IP: 1:500-1:5000

Specificity

Human, Mouse; other species not tested.

Antibody description

CD127/IL7R Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, K-562 cells. Positive IP detected in K-562 cells. Positive FC detected in K-562 cells. Observed molecular weight by Western-blot: 56-58 kDa

Preparation

This antibody was obtained by immunization of CD127/IL7R recombinant protein (Accession Number: NM_002185). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:109045(IL7R antibody) at dilution of 1:1500



IP Result of anti-IL7R (IP:Catalog No:109045, 4ug; Detection:Catalog No:109045 1:1000) with K-562 cells lysate 1200ug.



1X10⁶ K-562 cells were stained with 0.2ug IL7R antibody (Catalog No:109045, red) and control antibody (blue). Fixed with 90% MeOH blocked

CD127/IL7R antibody



Catalog Number: 109045

with 3% BSA (30 min). Alexa Fluor 488-

conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)
with dilution 1:1000.