

NUR77 antibody



Catalog Number: 113445

Product name

NUR77 antibody

Immunogen

[Human NR4A1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NUR77 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 70-80 kDa

Preparation

This antibody was obtained by immunization of NUR77 recombinant protein (Accession Number: NM_002135). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113445(NR4A1 antibody) at dilution of 1:100