

SFRP2 antibody



Catalog Number: 115145

Product name

SFRP2 antibody

Immunogen

[Human SFRP2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SFRP2 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, COLO 320 cells, MCF7 cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of SFRP2 recombinant protein (Accession Number: NM_003013). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:115145(SFRP2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115145(SFRP2 antibody) at dilution of 1:50 (under 10x lens)