

MEK4 antibody



Catalog Number: 112645

Product name

MEK4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MEK4 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 44-46 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003010). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112645(MAP2K4 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:112645(MAP2K4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:112645(MAP2K4 Antibody) at dilution of 1:50 (under 40x lens)