

# Anti-RTL6 antibody



Catalog Number: 176055

## Product name

Polyclonal

Anti-RTL6 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to RTL6

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ldoc1l. aa range:141-190

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1.; Western Blot analysis of various cells using LDOC1L Polyclonal Antibody