

PKIG antibody



Catalog Number: 113855

Product name

PKIG antibody

Immunogen

[Human PKIG Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

PKIG Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IF detected in Hela. Observed molecular weight by Western-blot: 5-8kd

Preparation

This antibody was obtained by immunization of PKIG recombinant protein (Accession Number: XM_047439849). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113855(PKIG antibody) at dilution of 1:200



Immunofluorescent analysis of Hela cells, using PKIG antibody Catalog No:113855 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).