

WDR13 antibody



Catalog Number: 116855

Product name

Polyclonal

WDR13 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

WDR13 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human heart tissue, Jurkat cells, mouse kidney tissue. Observed molecular weight by Western-blot: 43kd, 53kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of WDR13 recombinant protein (Accession Number: BC002507). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116855(WDR13 antibody) at dilution of 1:300

Clonality