

TNNT3 antibody



Catalog Number: 116155

Product name

TNNT3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TNNT3 Rabbit Polyclonal antibody. Positive WB detected in rat skeletal muscle tissue. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 35-40

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001042780). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116155(TNNT3 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:116155(TNNT3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:116155(TNNT3 Antibody) at dilution of 1:50 (under 40x lens)