

P5CS antibody



Catalog Number: 113555

Product name

P5CS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

P5CS Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, HepG2 cells, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 87 kDa

Preparation

This antibody was obtained by immunization of P5CS recombinant protein (Accession Number: NM_002860). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113555(ALDH18A1 antibody) at dilution of



IP Result of anti-ALDH18A1 (IP:Catalog No:113555, 3ug; Detection:Catalog No:113555 1:500) with mouse brain tissue lysate 5000ug.