

PLEKHA1 antibody



Catalog Number: 113955

Product name

PLEKHA1 antibody

Immunogen

[Human PLEKHA1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLEKHA1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human placenta tissue, Jurkat cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of PLEKHA1 recombinant protein (Accession Number: XM_017016481). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

Validations

230kd →
130kd →
95kd →
72kd →
56kd →
36kd →
28kd →
17kd →



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113955(PLEKHA1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113955(PLEKHA1 antibody) at dilution of 1:100 (under 25x lens)