

Kv1.2 antibody



Catalog Number: 112255

Product name

Kv1.2 antibody

Immunogen

[Human KCNA2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Kv1.2 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 66-75 kDa

Preparation

This antibody was obtained by immunization of Kv1.2 recombinant protein (Accession Number: NM_001204269). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations




WB result of anti-KCNA2 (Catalog No:112255) in



different cell and tissue lysates.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112255(KCNA2 antibody) at dilution of 1:1000



IP Result of anti-KCNA2 (IP:Catalog No:112255, 4ug; Detection:Catalog No:112255 1:500) with

Kv1.2 antibody



Catalog Number: 112255

HeLa cells lysate 5200ug.