

Anti-SF3B3 Rabbit antibody



Catalog Number: 167255

Product name

WB, ICC/IF, FC

Anti-SF3B3 Rabbit antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse, Rat

ICC/IF: 1/50

Antibody description

FC: 1/50

SF3B3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SF3B3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of SF3B3(green) in HeLa cells using SF3B3 Rabbit pAb at dilution 1/50, and



Western blot detection of SF3B3 in SHSY-5Y, C6 cell lysates using SF3B3 Rabbit pAb(1:1000 diluted). Predicted band size: 136KDa. Observed band size: 136KDa.